

## Activotec Peptide Synthesizer Users Give 5\*\*\*\*\* Reviews!

**This is a sample of what some of our customers have said about our synthesizers and support.**

### **University of Aalborg, Denmark**

'The Activo-P11 peptide synthesizer has been very successful in our laboratory producing quality peptides of good yields. The synthesizer has proved to be extremely reliable'

### **Tuebingen, Germany**

'I would like to say we are very satisfied with the Activo-P11 synthesizer, it was mainly the excellent service we experienced and the synthesizer reliability that led us to our decision to buy a second machine.'

### **IMM – Oxford, UK**

'The Activotec technical support has been excellent over the last 10 years, they have always fixed my peptide synthesizer quickly and the team are very knowledgeable about automated peptide synthesis.'

### **FMP – Berlin, Germany**

'the Activo P-11 has been a very reliable system, the possibility to run the entire synthesis or single steps at high temperature has been found to be very helpful for overcoming "difficult" steps. Very importantly, the personnel of Activotec are very supporting, and the service, if needed, is very quick and good.'

### **Cancer Research Institute, UK**

'An experienced PhD student has been using the Activo-P14 heavily since it was installed and has found it easy to use and versatile. His research has really moved forward as a result of being able to construct peptides at variable scales. I also have a new student started just 2 weeks ago, and already she's produced her first gram-scale peptide on the Activo-P14, so we're finding it good for what we want to do.'

### **Gdansk – Poland**

'We are really happy with the job that was done during the service visit. All is running good, and we are getting exceptionally pure peptides.'

### **Belgium**

'We purchased a P11 peptide synthesizer from Activotec in 2012. This equipment has been used successfully since then to synthesize peptides of varying length from 5 to around 80 amino acids. The synthesizer is equipped with UV monitoring and is easy to use.'

### **MPI – Germany**

'In Sept 2009, we bought an Activo-P11 peptide synthesizer including the UV monitor and the machine has been running almost constantly ever since. The set-up of the machine was very easy and the technical support for the installation of the machine was excellent. We are indeed very happy with the machine itself but also particularly with the services provided by the Activotec support team. Since we are not only doing classical peptide synthesis, we need to use unusual chemicals or change coupling protocols. In our experience it was relatively easy to do this with the Activo-P11 and so far we were able to perform all reactions on the automated system we wanted to. Besides we use the system also for classical peptide synthesis and have been particularly satisfied with the cheap and customer friendly supplies required such as the reactors and amino acid vials.'

### **Berlin – Germany**

'We were able to successfully synthesize even long and complicated sequences with the machine.'